

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



3128 B Blood Pressure Monitoring system

Measurement Range Systolic : 60-250 mmHg
Measurement Range Diastolic : 30-195 mmHg
Accuracy pressure : $\pm 2\%$

Pulse rate (PS) range : 0-199 slagen /min.

PS accuracy : $\pm 4\%$

Power source : 4 x AA batteries or charger

Weight : 350 gram

Data output : Bluetooth transmission or USB

For more information on this Blood Pressure Monitoring system please email to : info@prior-medical.com