

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



8201 Pulse Oxi Meter heavy duty

SpO2 display range	: 0 % tot 100 %
SpO2 resolution	: 1 %
SpO2 method	: dual wavelength LED
SpO2 Accuracy	: 80 tot 100%, ± 2 % 70 tot 79 %, ± 3 %
Heart rate (HR) range	: 30-250 beats /min.
HR accuracy	: ± 1 % or 1 beat/min if HR <100/min
Led Display	: SpO2, HR, 3 digits. Pulse strength bar
Alarm range setting	: high /low for SpO2 and pulse rate
Water resistance	: IPX2, watersplashes
Shockresistance	: 1 meter shockproof
Power source	: 4 x AA

For more information on this Pulse Oxi system please email to :
info@prior-medical.com