

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



8255 B Fingertip puls oximeter

SpO2 Range	: 0 % to 100 %
SpO2 Measurement methode	: dual wavelength LED
SpO2 Accuracy	: 80 to 100% , \pm 2 %
Perfusion index	: 5 bar indication
Heart Rate Range	: 30-250 beats /min.
Heart rate accuracy	: \pm 1 % or 1 beat/min
Communication via	: Bluetooth LE
Dinmensions	: 63 (L) x 35 (B) x 36 (H) mm
Fingertip surface	: silicon rubber
Power source	: 2 x AAA
Weight	: 40 gram
Auto off	: after 15 seconden

For more information email : info@prior-medical.com

