

Protocolbuster B.V.  
Ondernemingsweg 66m  
2404 HN Alphen aan den Rijn  
Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-085-0250122  
Email : [info@protocolbuster.com](mailto:info@protocolbuster.com)



**PROTOCOLBUSTER**



### **Bracelet tag, with alarm button.**

- |              |                                                                                                      |
|--------------|------------------------------------------------------------------------------------------------------|
| Weight       | : 23 grams                                                                                           |
| Class        | : IP 67                                                                                              |
| Material     | : Silica gel, anti-bacterial                                                                         |
| Colour       | : Grey                                                                                               |
| Shape        | : bracelet                                                                                           |
| Version      | : Bluetooth 4.2, compatible Bluetooth 5                                                              |
| Range        | : 40 meters                                                                                          |
| Battery life | : 9 months at 1000 ms sending (12 hours powersafe)<br>5 months at 350 ms sending (12 hour powersafe) |
| Battery type | : CR 2032 Lithium Manganese Dioxide Battery                                                          |

For more information: [info@prior-medical.com](mailto:info@prior-medical.com)  
This product is delivered in cooperation with Protocol Buster. This product is used in combination with the location manager.