

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



PMG 80, digital dynamometer for physical therapy and recovery.
Including case.

- Range : 20/ 40 / 80 kg
- Graduation : 100 gram
- Power supply: button battery 3v type CR2450 li
- Functions : counting, average right/left, maximum
- Grip size : 64 mm

For information: info@prior-medical.com