


Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model PMS 2510, Scale for lifters to weigh immobile patients.

- Capacity **300 kg**
- Graduation per **100 gram**
- **Functions : hold en auto off**
- Digital scale by loadcell
- Dimensions : 52 x 122 x 160 mm
- Power : battery
- **OIML, Class III calibrated, MDD 93/42/CEE amendment 2007/47/CE**

For information : info@prior-medical.com