

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Model PMS 3830 , Multifunctional Wheelchair and Stretcher Weighing Scale
for hospitals and clinics, dialyses centre and long stay care.

- Weighing capacity **300 kg**
- Graduation per **100 gram**
- Functions : BMI, Hold, Tare
- RS 232 Port for printer /computer connection
- Fully electronic with 4 OIML approved loadcells
- Movable through 2 wheels at the side
- Weight : 34,5 kg
- Platform drive in dimensions 900 (l) x 754 (b) mm
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**

For more information email : info@prior-medical.com

