

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model PMS 4400, Handheld baby scale. For mobile situations.

- Graduation per 10 gr
- Functions; on/off, tare (100% of capacity) and hold (if baby moves a lot)
- Maximum capacity **10 kg**
- Clear display (25 mm)
- Including 6 batteries
- Dimensions 79x150x193 mm
- **Calibrated, Class III, MDD 93/42/CEE amendment 2007/47/CE**

For information : info@prior-medical.com

