

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



**Prior Model PMS 5461**, Movable Multifunctional Chair Weighing Scale for weighing patients

- Weighing capacity **250 kg**,
- Graduation per **100 gram**
- Fully electronic loadcell weighing
- With foldable armrests
- For the stability during weighing of the patient equipped with 2 foldable feet
- Chair standard movable by 4 castors, with a brake on it
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**

For information : [info@prior-medical.com](mailto:info@prior-medical.com)