

Prior Medical Systems B.V.

Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086

Email : info@prior-medical.com



Prior Medical Systems



PMS 5710

Multifunctional Electronic Weighing Scale

- Maximum capacity 300 kg
- Graduation per 100 gram
- Functions; on/off, tare (100% of capacity)
- Clear display (29 mm)
- Automatic on/off
- Including adapter and 6 batteries
- Dimensions Platform 310 x 360 mm , including 4 adjustable feet and level indicator
- Handle
- **Calibrated, Class III, MDD 93/42/CE amendment 2007/47/CE**

For information : info@prior-medical.com

