

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



**Prior Model PMS 5900**, Multifunctional Baby scale to weigh quickly, easy and safely baby's and infants

- Graduation per 2 **gram** (0 ~ 6 kg) and 10 **gram** (6 ~ 15 kg)
- Functions; on/off, tare (100% of capacity) and hold (if baby moves a lot)
- Maximum capacity **15 kg**
- Clear display (20 mm)
- Automatic on/off
- Including adaptor and 6 batteries
- Dimensions 660 x 334 x 125 mm , including 4 adjustable feet and level indicator
- Optional digital rod
- **Calibrated, Class III, MDD 93/42/CEE amendment 2007/47/CE**

For information : [info@prior-medical.com](mailto:info@prior-medical.com)