

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



**PTD 1090**, for the calibration of the TD 1261 and 1107 thermometers, automatically recalibrates TD 1107. Includes an aluminium carrying case.

Verification of calibration accuracy : TD 1107 and TD 1261  
Auto recalibration : TD 1107

Power and fuse : 100-240 V (50.60 Hz/ 3  
Operating conditions : 16-32 °C.  $\leq$  85% RH  
Memory : 1000  
Temperature range : 37.0 – 41.0 °C  
Adjustable recalibration range :  $\pm$  0,5 °C  
Warm up time : 30 minutes  
Computer interface : USB HID cable

Dimensions : 300 x 350 x 160 mm  
Weight : 4,3 kg

For information or a quote, contact via email : [info@prior-medical.com](mailto:info@prior-medical.com).