

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model P BW 0365, Multifunctional Baby scale to weigh quickly, easy and safely baby's and infants

- Graduation per 2 **gram** (0 ~ 6 kg) and 10 **gram** (6 ~ 15 kg)
- Functions; on/off, tare (100% of capacity) and hold (if baby is moving)
- Maximum capacity **15 kg**
- Clear display (20 mm)
- Automatic on/off
- Including adapter
- Dimensions 600 x 300 x 160 mm , including 4 adjustable feet and level indicator
- Optional mechanical height rod, RS 232 output
- **Calibrated, Class III, MDD 93/42/CEE amendment 2007/47/CE**

For information : info@prior-medical.com