


Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



## **Prior Model P BW 0520, Multifunctional Electronic Weighing Scale**

- Maximum capacity 200 kg
- Graduation per 100 gram
- Functions; on/off, 0
- Removable and mountable display
- Platform with grip
- Clear LCD display
- Automatic on/off
- Including adapter
- Dimensions 340 x 330 x 66 mm , including 4 adjustable feet and level indicator
- **Calibrated, Class III, MDD 93/42/CE amendment 2007/47/CE**

For information : [info@prior-medical.com](mailto:info@prior-medical.com)

