

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Fax: +31-172-766105
Email : info@prior-medical.com



Prior Medical Systems



Model Prior 2000 / BW 1430 Multifunctional Electronic Weighing Scale with handrail for unstable or obese patients. With mechanical rod.

- Maximum capacity 300 kg
- Graduation per 100 gram
- Functions; on/off, auto off, hold, tarra and BMI
- Dimensions platform: 600 mm (l) x 600 (w) mm
- Weight: 19 kg
- Automatic on/off
- Including adaptor and batteries
- Movable through 2 wheels at the side
- With handrail and a mechanical length meter 110-200 cm / 0.1 cm
- **Calibrated, Class III**

For information : info@prior-medical.com