


Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Model Prior 3001, Movable Multifunctional Chair Weighing Scale for the weighing of patients

- Weighing capacity **300 kg**,
- Graduation per **100 gram**
- Fully electronic loadcell weighing
- With foldable armrests
- Dimensions seating : 48 x 50 cm
- For the stability during weighing of the patient equipped with 2 foldable feet
- Chair standard movable by 4 castors, of which 2 have a brake
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**

For information : info@prior-medical.com