

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model 3600, Multifunctional **Electronic Weighing Scale with integrated height rod and automatic BMI calculation** .

- Maximum capacity **250 kg**
- Graduation per **100 gram**
- Functions; on/off, tare (100% of capacity) and automatic BMI
- Integrated height rod 110-200 cm
- Clear display (29 mm)
- Automatic on/off
- Including adaptor
- Dimensions Platform 320 x 375 mm , including 4 adjustable feet and level indicator
- Option : height rod
- **Calibrated, Class III, MDD 93/42/CE amendment 2007/47/CE**

For information : info@prior-medical.com