


Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior 4000 W

Multifunctional **Electronic Weighing Column Scale** .

- Maximum capacity 250 **kg**
- Graduation per 100 **gram**
- Functions: Hold, BMI
- Clear display (29 mm)
- Automatic on/off
- Including 6 batteries and adaptor
- Dimensions 34.0 x 34.0 x 15,7 cm
- Including 4 adjustable feet and level indicator
- Wifi transmission
- **Calibrated, Class III**

For information : info@prior-medical.com