

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : info@prior-medical.com



# Prior Medical Systems



**Model Prior 4500** , Multifunctional Wheelchair Scale for hospitals and clinics, dialysis centres and long stay care.

- Weighing capacity **300 kg**
- Graduation per **100 gram**
- Functions : BMI, Hold, Tare
- Powered by battery or adapter
- Fully electronic with 4 OIML approved loadcells
- Movable through 2 wheels at the side
- Weight : 42 kg
- Platform drive in dimensions 760 (l) x 800 (w) mm
- Foldable handle with mounted indicator
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**