

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Model Prior 7000 , Multifunctional Wheelchair Scale for hospitals and clinics, dialysis centres and long stay care.

- Weighing capacity **400 kg**
- Graduation per **200 gram**
- Functions : BMI, Hold, Tare
- Powered by battery or adapter
- Fully electronic with 4 OIML approved loadcells
- Movable through 2 wheels at the side
- Weight : 52 kg
- Platform drive in dimensions 1000 (l) x 1000 (w) mm
- Handgrips 3 sides
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**