

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior BW 6190 class III wheelchair/bed scale for weighing non-mobile patients.

- Weighing capacity up to 360 kg
- Weight division per 100 grams
- Entry and exit
- Weighing section dimensions 1000 (w) x 1900 (l) mm
- Total entrance and exit dimensions: 1270 (W) x 2380 (l) x 112 (h) mm
- Own weight approx. 111,2 kg
- Calibrated, Class III, MDD 93/42/EEC Amendment 2007/47/EC

If you have any questions about this scale Prior BW 6190, please do not hesitate to contact us by email : info@prior-medical.com