

Prior Medical Systems B.V.

Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-84-4300044
Fax: +31-172-766105
Email : info@prior-medical.com



Prior Medical Systems



Prior PW 1536, trolley of synthetic material for baby scales. Worktop with raised edge. Drawers can be pulled out 110% by telescopic guides. Noiseless wheels.

- Maximum load capacity: 55 kg
- 4 braked wheels Ø 100 mm
- Push handle on 2 sides
- synthetic frame
- Color drawers: blue
- Drawers: 2, drawer load 15 kg max
- Dimensions: External 970 x 430 x 840 MM (l x w x h)
Work surface 660 x 410 MM
Drawer 560 x 360 x 120 MM

For more information : info@prior-medical.com