

Prior Medical Systems B.V.

Ondernemingsweg 66M
2404 HN Alphen aan den Rijn

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior PW 1773, trolley for baby scales. Worktop with raised edge. Drawers can be pulled out 110% by telescopic guides. Noiseless wheels.

- Maximum load capacity: 75 kg
- 4 braked wheels Ø 100 mm
- Push handle on 1 side
- Stainless steel frame
- Color options for drawers: green, blue or red
- Drawers: 3, drawer load 25 kg max
- Dimensions: External 730 x 420 x 1030 MM
Work surface 680 x 420 MM
Small drawer 560 x 360 x 120 MM
Large drawer 560 x 360 x 240 MM

For more information: info@prior-medical.com