

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Ear infrared Thermometer TD 1107, fast reaction time (2 seconds), compact

Temperature range	: 32° C – 43° C
Method	: ear
Accuracy	: ± 0.2 °C
LCD back light	: colours green when normal
Probe cover detection	: activates auto on function
Communication	: bluetooth
Power source	: 2 x AA batteries

Automatic calibration with the TD 1090 Calibrator

For more information on this thermometer please email to :
info@prior-medical.com