

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Thermometer TD 1241 B, infrared, fast reaction time, forehead measurement

Measuring distance : 3 – 7 cm
Range : 22 °C – 44 °C
Measurement method : forehead, infrared focus
Accuracy : 36 °C – 39 °C \pm 0.2 °C. (< 36 °C > 39 °C \pm 0,3 °C)
Resolution : 0,1 °C
Communication : bluetooth
Power source : 2 x AA batteries
Dimensions : 150 mm(l) x 48.5mm (b) x 55mm(d)
Weight : 125 gr (incl batteries)

Conformiteit : 93/42/EEC, EN 60101, EN 61326,
EN 62304

Clinically tested, results available on request.

For more information on this thermometer please email to :
info@prior-medical.com