

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Thermometer TD 1242 H, infrared, fast reaction time (1 second), forehead measurement

Measuring distance	: 2 – 5 cm
Range	: 32 °C – 43 °C
Measurement method	: forehead, infrared focus
Measurement modes	: baby, adult, surface
Display	: flashes with fever
Accuracy	: ± 0.2 °C, 35,0°C – 42,0 °C $\pm 0,3$ °C, < 35,0°C -> 42 °C
Resolution	: 0,1 °C
Power source	: 2 x AAA batteries
Dimensions	: 155.46mm x 40.14mm x 39.45m
Weight	: 61,8 gr (without batteries)
Auto off	: 30 seconds
Conformity	: 93/42/EEC, EN 60101

For more information on this thermometer please email to :
info@prior-medical.com