

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



## Ear infrared Thermometer TD 1261, fast reaction time (2 seconds), connects with TD 2300 monitor

Temperature range	: 32° C – 43° C
Method	: ear
Accuracy	: $\pm 0.2$ °C
Communication	: IrDA communication with 2300
Probe cover detection	: only works with probe cover
Power source	: 2 x AAA batteries
Dimensions	: 160.5mm(b) x 45.9mm(b) x 32.5mm(d)
Weight	: 98, 33 gr

For more information on this thermometer please email to :  
[info@prior-medical.com](mailto:info@prior-medical.com)